



ENHANCED WOUND HEALING
WITH **POLYMEM**® SILICONE BORDER

12 OCTOBER
2021

SHARON GARDNER
PROF. LINDA RAFTER

LUNCHTIME
13:00 - 13:45

Learning outcomes

1

Develop understanding of current wound healing challenges.

2

Recognise what makes **PolyMem** unique.

3

Learn how **PolyMem** can benefit your patients.

4

Learn about the **PolyMem** dressing range.

Wound healing challenges


- Clinicians are under pressure to provide high quality patient outcomes at reduced cost
- Cases presented are complex wounds that require careful management
- Some patients were in there 90s, extremely compromised with nutritional deficiencies
- Honey and superabsorbent dressings were used prior to the patients being commenced on the **PolyMem** dressings.

What is PolyMem?

PolyMem is a unique multifunctional polymeric membrane dressing.

It consists of:

- A mild, non-toxic wound cleanser (surfactant)
- A soothing moisturiser (glycerol)
- A superabsorbent (starch co-polymer)
- A semi-permeable film backing.*




*Not included in WIC dressings.

How does it work?



How does it work?



Key benefits


Modulates
inflammatory response

Relieves pain

Supports wound
healing



Structure of nociceptors



Key benefits

- Wound cleanser is released continually into the wound, loosening bonds between wound debris and healthy tissue
- The moisturiser (glycerol) is simultaneously released, creating a moist wound environment, preventing the dressing from sticking
- **PolyMem** supports autolytic debridement.



Case study one: debridement





- 63-year-old female
- Anaesthetic risk due to past medical history
- Sustained large haematoma tibia/calf area
- Wound bed preparation with **PolyMem**
- Larval therapy
- **PolyMem** to wound healing.

Indications

- Abrasions/trauma
- Bruising
- Superficial burns
- Skin tears
- Surgical wounds
- Diabetic foot ulcers
- Fungating wounds
- Leg ulcers
- Pressure ulcers
- Dermatology e.g. *epidermolysis bullosa* (EB)
- Donor and graft sites
- Exposed tendons
- Radiotherapy skin reactions.

PolyMem MAX for Pronton therapy



How do we use it?

- Dressing activation
- Cutting guide – **coming soon!**



Application video

Video Here

Case study two: trauma wound



- 73-year-old female
- Cardiac and oncology history
- lady fell and sustained a trauma wound 05/12/19
- this led to the right tibia with significant tissue loss and bone exposure
- 03/02/20 she commenced treatment with **PolyMem** dressings.

Case study two: outcome



- Final assessment on 06/03/20, wound was healed
- Using **PolyMem** dressings for **32 days** accelerated this lady's healing and prevented her requiring a skin graft.

NEW PolyMem Silicone Border

- The same patented polymeric membrane formulation with a secure, gentle adhesive
- Soft silicone border is gentle to the skin and repositionable
- Designed to **facilitate healing, relieve pain and reduce inflammation**
- An ideal choice for patients with fragile, sensitive skin.



Case study three: skin tear



- 14/06/21, female sustained a Skin Tear Audit Research (STAR) classification 2b (Carville et al, 2007)
- There was slight swelling to the left lower leg
- Pain score on the McCaffery scale (Pasero and McCaffery, 2011) was 9
- 05/07/2021 commenced treatment with **PolyMem** dressings.

Case study three: outcome



After 7 days using **PolyMem MAX**:

- Purple tissue had resolved
- Pain level had reduced to 0
- The wound had reduced significantly in size.

NEW PolyMem Silicone Border

- Suitable for dry-to-moderate exuding wounds
- Suitable for very fragile skin
- Easy to apply
- Very easy to remove without any trauma to the skin.



Case study four: trauma



12/07/2021 **PolyMem Silicone Border** was applied


- Wound had reduced in size by 1x0.5cm
- Exudate volume reduced
- The dressing was easy to apply
- Dressing changed every seven days
- The wound healed in five weeks.

Case study five: pressure ulcer



- 88-year-old female
- Presented with multiple pressure ulcers
- Significant past medical history, including dementia
- Poor oral intake
- Challenges with repositioning
- February 2021 thoracic spine pressure ulcer commenced on **PolyMem**.

Case study five: outcome



28/09/2021 — pressure ulcer has nearly healed

- **PolyMem Silicone Border** was applied
- Stayed in place and easy to remove without causing any damage to very fragile skin
- Her pain level was assessed by facial expressions due to dementia.

PolyMem Finger dressing



- 39 patients were evaluated
- Standard dressing was compared to the **PolyMem Finger** dressing for wound healing by secondary intention in trauma injuries
- Treated in general hospital accident and emergency and minor injuries units.

PolyMem Finger dressing: aims

- Explore the patients' perceptions of pain by using numerical and descriptive pain diaries
- Explore how patients' quality of life was affected: hygiene needs, dressing, ability to work, driving, domestic chores, mobility and hand function
- Monitor the patients' sleep patterns in line with injury and pain
- Explore patients' pain at dressing change
- To compare the costs and procedure time using the two dressings evaluated
- Monitor the patient use of over-the-counter analgesic
- To explore the nurses' perceptions of the dressing change.

PolyMem Finger dressing




- 24-year-old male
- 07/11/13 — trapped finger in car door
- Surgical debridement
- **PolyMem Finger** dressing
- 25/11/13 — healed with minimal scarring
- **Pain score reduced from 10 to 0 within two weeks of injury.**

PolyMem: the range

- NEW PolyMem Silicone Border
- PolyMem Film Island
- PolyMem Non adhesive
- PolyMem MAX
- PolyMem MAX Film Island
- PolyMem Finger/Toe
- PolyMem Tube
- PolyMem Silver.



PolyMem MAX



- For greater absorbency and longer wear
- This product is 60% thicker than original **PolyMem**
- It is a specialised dressing with multiple properties:
 - Moisturises
 - Cleanse and debrides
 - Analgesic properties
 - More absorbent than standard **PolyMem**.

Summary of findings

- **PolyMem MAX** can be used to debride a wound bed and relieve pain
- **PolyMem Silicone Border** is extremely gentle on the skin and stays in place well
- Removed very easily without trauma
- All patients had particularly good outcomes
- Case studies included very complex wounds for which it is particularly difficult to gain complete healing.

Conclusion


- Clinicians need to work together with their patients to promote cost-effective wound healing
- All case study patients were involved in their care pathways, and it was explained how **PolyMem** promoted wound healing.

Ensuring the **right dressing** is used at the **right time** on the **right patient** can ultimately result in **good patient outcomes**.

Educational support

Facebook live with GPN

Understanding wound pain: the physiology of nociceptors and inflammation.



E-learning module with Wound Care Today

Identifying and managing wound pain.



Visit our 'Education' tab here: www.hrhealthcare.co.uk

Clinical support

- PolyMem patient and user guides
- YouTube application videos
- NEW published paper British Journal of Community Nursing 'Achieving effective patient outcomes with *PolyMem Silicone Border*', by Linda Rafter, Mark Rafter
- Cutting guide – coming soon!
- PolyMem brochure
- PolyMem Silicone Border brochure.



Get in touch!

If you would like to evaluate **PolyMem Silicone Border** in your organisation,
please email:



Marketing@hrhealthcare.co.uk

Or, get in touch with your local representative for more information.

References

- Carville K, Lewin G, Newall N, et al (2007) STAR: a consensus for skin tear classification. *Prim Intention* **15(1)**: 18-28
- Pasero C, McCaffery M (2011) *Pain Assessment and Pharmacologic Management*. Mosby-Elsevier, St. Louis
- Rafter L, Oforka E (2013) Trauma-free fingertip dressing changes. *Wounds UK* **9(1)**: 96-100
- Rafter L, Rafter M. (2021) Achieving effective patient outcomes with PolyMem dressing. *Br J Community Nurs* **26(10)**: 2-10

Discussing: **Polymem**® Silicone Border





Download your certificate

www.jcn-live.co.uk/virtual-lunch-meeting